

• 临床检验研究论著 •

Validation and quality control of hematolymphoid neoplasm immunophenotyping by flow cytometry

Wu Lijuan¹, Xu Dongsheng²

(1. Department of Clinical Laboratory Medicine, Chengdu Military General Hospital, Chengdu, Sichuan 610083, China; 2. Department of Hematopathology, CBLPATH, Inc. Rye Brook, New York 10573, USA)

Abstract: cytometric immunophenotyping has evolved from two-parameter quantitative measurement of peripheral blood lymphocytes to five-or more parameter qualitative evaluation of bone marrow for hematopathology. Leukemia/lymphoma immunophenotyping represent an important addition to histomorphology in the diagnosis, classification and monitoring of hematolymphoid neoplasms. The complexity of five- or more parameter analyses and the interpretation of the data rely on standardization and validation of the instrument, the reagent and the procedure. In addition, clinical flow cytometry laboratories in U. S. are required to document proficiency testing, sample preparation, method accuracy, specificity, sensitivity and precision. CLSI and the U. S. -Canadian Consensus Conference have provided recommendations, but each laboratory is responsible for validating its own qualitative and quantitative procedures. This paper introduces the procedures for quality control of all levels of the operation in a clinical flow cytometry laboratory in USA.

Key words: flow cytometry; immunophenotyping; quality control

DOI: 10. 3969/j. issn. 1673-4130. 2012. 23. 018

文献标识码: A

文章编号: 1673-4130(2012)23-2852-07

Clinical flow cytometry has evolved from two-parameter quantitative assessment of peripheral blood lymphocytes to five-or more parameter qualitative evaluation of bone marrow for hematopathology. Leukemia and lymphoma immunophenotyping represent an extremely important complement to morphology in the diagnosis, classification and monitoring of hematolymphoid neoplasms.

The complexity of five- or more parameter analyses and the interpretation of the data rely on standardization and validation of the instrument, the reagents and the procedure. In addition, flow cytometry laboratories in the U. S. are required to document proficiency testing, sample preparation, method accuracy, specificity, sensitivity and precision. CLSI and the U. S. -Canadian Consensus Conference have provided recommendations, but each laboratory is ultimately responsible for validating its own qualitative and quantitative procedures. This paper reviews procedures for validation and quality control of all aspects of the operation of a clinical flow cytometry service.

1 Introduction

Flow cytometry is a dynamic technology which has allowed the multi-parameter analysis of heterogeneous cell populations to develop as a clinical service. Complex analyses are able to combine immunophenotyping of both surface and cytoplasmic antigens, DNA analysis and functional evaluations. Subsets of cells can be identified and characterized by patterns of maturation antigens and staining intensity which can assist in diagnostic and prognostic interpretations as well as the detection of minimal residual disease^[1].

Four- or more color immunophenotyping are not unusual in the clinical flow cytometry laboratory these days, allowing the simultaneous measurement of multiple different parameters. With many variables in these analyses, standardization and validation of instrumentation and methodology is essential to

ensure the technical quality of the results. Reagents must be well characterized for specificity and performance with different fluorochromes. All monoclonal-fluorochrome combinations must be critically evaluated for staining intensity, spectral overlap, and instrument compensation. Furthermore, fluorescence patterns must be characterized for diagnostic combinations of antigens and diagnostic interpretation.

There is growing recognition of the need to set standards of training and education for practitioners in the field of flow cytometry, including technologists, interpreters, and laboratory directors. In attempts to assist with standardization, the U. S. -Canadian Consensus Conference in 1997^[2] and Clinical and Laboratory Standards Institute (CLSI) guidelines in 2007^[3] provided recommendations for clinical flow cytometry in hematopathology. On the regulatory side, the Clinical Laboratory Improvement Act (CLIA'88) has influenced laboratory staffing, training, validation and documentation. As we enter the 21st century, the laboratory's responsibilities continue to increase.

The goal of this article is to discuss how regulatory oversight influences the laboratory's validation and quality control documentation, particularly in hematopathology. Although the article covers regulatory issues applying in the U. S., all clinical laboratories world-wide are responsible for maintaining performance standards. Good Laboratory Practices include Standard Operating Procedures (SOPs), Quality Control (QC) and Quality Assurance (QA) as integral to providing good patient care. Regardless of requirements placed on manufacturers to provide analyte-specific reagents (ASR) or in vitro diagnostic (IVD) reagents, each laboratory should validate their own panels for sensitivity, specificity, and correlation with morphology and clinical findings. The following sections will consider laboratory responsibilities, reagent and instrument validations, QC, QA, and troubleshooting as a model for the integration of new

technologies into diagnostic hematopathology.

2 Responsibilities of the clinical flow cytometry service

Flow cytometry falls under the CLIA'88 category of high complexity laboratory testing. Documentation of staff qualifications and training as well as analytical accuracy, sensitivity, precision and QC are required. In some states, such as California, personnel responsible for generating flow cytometric results must be licensed medical technologists who are required to have 12 h of continuing education yearly. In New York, the flow results must be reviewed by an individual holding a certificate of qualification issued by New York State Department of Health^[4]. At the very least, training and proficiency in the technology must be documented. Laboratories should develop SOPs for training, with supervisor verification of proficiencies in incremental responsibilities for each staff member. Instrument manufacturers and scientific societies also conduct training courses and provide certificates of training. Additionally, the American Society of Clinical Pathology (ASCP) offers a specialty exam for national certification in flow cytometry.

2.1 Proficiency testing All high-complexity laboratories must enroll in a proficiency testing (PT) program that meets CLIA regulations and is approved by Department of Health and Human Services (HHS). Regulations require that the laboratory inform HHS which programs it will use, list the tests performed for each program, participate in a program for at least a year before choosing a different program, and must notify HHS before changing. The laboratory must also authorize the PT program to release data to HHS.

Proficiency Testing samples must be treated the same as patient specimens, using the same personnel and work processes. Results must not be discussed with other laboratories and PT samples may not be outsourced under any circumstances. Documentation for handling, preparation, processing, analysis and interpretation of PT samples must be kept. Final reports must be retained for a minimum of 2 years^[4].

2.2 Accuracy Good laboratory practices and CLIA'88 require that all clinical testing be characterized for accuracy, specificity, sensitivity and precision. However, for leukemia and lymphoma immunophenotyping by flow cytometry, there are no consensus standards or recommendations for these assessments. Thus each laboratory must establish its own performance criteria.

Analytical accuracy compares the test result to a reference, or 'gold' standard. In hematology, normal cell populations can be counted microscopically with a hemocytometer or with automated hematology analyzer and results compared for accuracy. In abnormal populations, especially when the characterization is based on multi-parameter data, the comparison becomes complex. For hematopathology, the gold standard is morphology. The flow cytometric assessment of accuracy must therefore compare to morphology. With rare events, however, morphologic diagnosis is difficult. Flow cytometry allows characterization of cell populations with complex phenotypes and should also assess accuracy by comparison to previously char-

acterized cells. Sources of well-characterized cell populations are cryopreserved pedigreed specimens from another validated laboratory. Another source of specimens for assessing accuracy can be cases diagnosed by cytogenetics and molecular biology. Documentation comparing at least 10 leukemia and lymphoma cases with complete histopathology and clinical findings should be on file to support analytical accuracy^[2,5].

Another aspect of accuracy that should not be ignored is the qualitative staining pattern used to identify cell lineages. Descriptions of dim, moderate, and bright staining patterns can be diagnostic or prognostic and should be well characterized and documented. For example, when multiple antibodies conjugated to fluorochromes are used in combination, accuracy should be assessed for each antibody separately and the frequency and intensity results compared to those obtained when using the combined reagents.

2.3 Specificity Specificity of monoclonal reagents is defined by how well the antibody recognizes the correct antigenic target. Manufacturers are responsible for reagents having the correct specificity listed on their labels. However, specificity in a 'home-brew' diagnostic test has a broader meaning. For leukemia and lymphoma testing, the interpretation of the flow panel should be compared with morphology and clinical presentation to assess the 'specificity' of the testing. Each laboratory should establish its own expected rate of discrepancy between flow and morphology, most likely, 5%. The laboratory should then assess, on a case-by-case basis, the reason for a discrepancy, document the discrepancy as a QA assessment, and monitor trends quarterly^[2-3,6].

Specificity of flow cytometric reagents can be assessed by consensus workshops, publications or by in-house validation. Leukocyte Typing VI aptly and succinctly summarize the vast body of testing results on most of the monoclonal antibodies commonly used in leukemia and lymphoma immunophenotyping. Other sources of specificity data are clinical texts and journal articles or email communications through Purdue Cytometry Mailing List.

However, flow cytometry monoclonal reagents are considered 'home-brew' reagents and each laboratory must document staining performance, or 'diagnostic specificity' on both normal and abnormal cell populations, identifying both positive and negative staining. Clinical hematology specimens will almost always contain normal along with any abnormal populations and thus will allow both evaluations in one specimen. Diagnostic immunophenotyping should be validated with pedigreed malignant specimens or well-characterized cell lines used as positive controls.

2.4 Sensitivity Reagent sensitivity describes the ability to detect the minimum staining intensity above non-specific or negative staining. The sensitivity of detection is dependent on the titration of monoclonal reagents, the proper instrument setup and calibration, the number of cells counted and the flow rate of the instrument. Documentation of reagent titrations and parallel testing of each new lot of antibody is required. Instru-

ment calibration and documentation is also required.

2.5 Precision Precision is a standard analytical parameter which measures the reproducibility of a single sample stained and analyzed in duplicate at least 10 times. Normal peripheral blood, cell lines or CD Chex may be used. CD Chex is preserved white cell controls manufactured by Streck Laboratories. Quantitative mean and standard deviation for each monoclonal antibody should be documented and a 2 S. D. range developed. Determination of instrument precision is accomplished by running the same stained sample at least three times, with results within 2 S. D.

2.6 Analyte-specific reagents In the United States, clinical laboratories, regulated under CLIA '88 and performing physician ordered flow cytometric testing for hematolymphoid immunophenotyping, are required to be qualified to perform high-complexity testing and must use reagents that are manufactured and labeled as Analyte Specific Reagents (ASRs).

Manufacturers of ASR must register their facility, list their products, follow current Good Manufacturing Practices (cGMP) under the new Quality System Regulations (QSRs), and comply with medical device reporting requirements. Manufacturers are restricted from providing any statement regarding analytical or clinical performance, number of tests provided, or instructions for use. Advertising and promotional material for an ASR product must include the identity and purity (including source and method of acquisition) of the analyte-specific reagent and the identity of the analyte.

The identity of the analyte for monoclonal antibody reagents is provided by the antigen distribution and supporting references. The labeling for analyte-specific reagents must include the reagent name, concentration, purity and quality, statement of warnings or precautions for users, date of manufacture, lot number, expiration date, storage instructions, net quantity of contents, name and place of business, and the following statement: "Analyte Specific Reagent, Analytical and performance characteristics are not established".

It is the responsibility of the testing laboratory using these reagents to validate their performance in 'home-brew' clinical assays. In addition, any laboratory reports using ASRs must contain the following disclaimer: "This test was developed and its performance characteristics determined by (laboratory Name). It has not been cleared or approved by the U. S. Food & Drug Administration"^[3-4]. Further clarifying comments may also be included.

Research institutions and other non-clinical laboratories that use ASRs to make tests for purposes other than providing diagnostic information to patients and practitioner are not restricted under the ASR regulation.

3 Validation of Immunophenotyping

The basic components are the validation of the instrument, the individual reagents, and the staining protocols used to create a final interpretation. Each step in the testing process must be specified and quality control measurements included monitoring performance.

3.1 Instrument validation Validation of instrument performance falls into two areas: (a) instrument setup and daily qualification of both light scatter and fluorescence measurements, and (b) cross-instrument performance using relevant clinical specimens.

3.2 Instrument setup and daily qualification CLSI recommends that the setup of a flow cytometer is comprehensive enough to assure proper optical alignment for adequate light scattering and fluorescence sensitivity and resolution as well as the proper degree of compensation to correct for spectral overlap when multiple fluorochromes are used. While hardware configuration may not allow the laboratory to align optics, a check of performance must be made and documented. In addition, since the hardware and optics are different for flow cytometers made by different manufacturers (the two major ones being Beckman Coulter and BD Biosciences), the recommended procedures for instrument setup and performance qualification are instrument-specific^[3-4].

3.3 Light scatter sensitivity and resolution Optical alignment for optimal sensitivity and resolution of both forward (FSC) and side (SSC) scatter can be assessed by running uniformly sized beads that fall within the light scatter ranges observed with most clinical samples, at a constant photomultiplier tubes (PMT) voltage on a daily basis. The mean FSC and SSC channel numbers and percent coefficient of variation (% C. V.) should be recorded. The acceptable ranges for each parameter can be established by first running the beads 20 times over a 5 day time period at the same PMT setting. Levy-Jennings graphs are then used to plot the values obtained daily and an action plan is established for what to do when any parameter falls outside of the expected range. However, since beads and cells do not always behave similarly on a flow cytometer, it is also recommended that the instrument be set up for running clinical samples using biological material and that the proper resolution of different cell types be determined in daily instrument qualification procedures. Normal peripheral blood leukocytes processed in a similar manner as the clinical samples can be used for this purpose^[5-6].

3.4 Fluorescence sensitivity and resolution Acceptance values for monitoring fluorescence sensitivity and resolution can be established by two methods: (a) recording the channel number and C. V. of calibration beads with a pre-determined, fixed laser power, filters, PMT voltages and gains, or (b) recording the high voltage and gain to position the beads in the same channel each time. For either approach, acceptance ranges can be established by using at least 20 replicate data sets collected over at least a 5 day period; new ranges must be developed each time a new lot of beads are used. Regardless of the method applied, procedures for what to do if any of the measured and plotted parameters do not fall within the acceptable range must also be established, followed and documented by the clinical laboratory. Values for the acceptable ranges displayed on the Levy-Jennings graphs for each parameter are typically updated based on each successive 20 day cycle of data collection. However, QC

specifications for instrument performance must be established to allow for detection of significant trends or drifts over time which require corrective action^[6].

3.5 Fluorescence compensation Compensation should also be evaluated at the time of initial instrument setup and then monitored daily, preferably using a biological compensation control in addition to hard dyed beads provided by the instrument manufacturer. Depending on the number, type and combination of fluorochromes used in a particular staining protocol, the compensation between fluorescence signals must be properly set to prevent spill-over of one fluorescence signal into another. This is usually done automatically using hard dyed beads for four-color compensation of FITC, PE, PerCP and APC on a FACSCalibur (BD Biosciences).

For six or more color flow cytometer, such as Gallios (Beckman Coulter), daily instrument optimization procedures include running Flow-Check Pro Fluorospheres (a suspension of fluorescent microspheres used for the daily verification of the cytometer optical alignment and fluidics system), the Flow-Set Pro Fluorospheres (a suspension of fluorescent microspheres used in the standardization of light scatter, fluorescent intensity, and hydrodynamic focusing to ensure optimal instrument performance on a daily basis), and the Cyto-Comp cell kit (consisting of lyophilized human lymphocytes and used to adjust color compensation settings prior to multi-color analysis. The cyto-comp cells with CD45 in every fluorochrome establishes a compensation matrix and the verify tube (CD8 FITC/CD4PE/CD3ECD/CD56PC5. 5/CD7PC7/CD20APC/CD19 APC-Alexa Fluor 750/ CD 45 Pacific Blue) as well as confirms the instrument settings. Once these beads have been acquired, the results are approved and stored in a Levy-Jennings file which can be printed on a monthly basis and filed in a QC binder.

4 Monitoring instrument performance

In the setup and qualification of a flow cytometer for daily use, CLSI recommends a check of optical alignment, fluorescence resolution and intensity. For the Coulter Gallios, Flow-Check Pro Fluorospheres is run daily to verify fluidics system and a proper optical alignment between the laser beam and the cell at the interrogation point as well as adjustment of laser power and/or PMT voltage.

For the BD FACS Calibur, the fluidics and calibration check is performed by running CaliBrite Beads and allowing the instrument software to make adjustments in setup based on lot-specific performance expectations. Documentation of acceptable calibration is provided for record keeping. These records can be used for trend analysis of the laser power and voltage settings to indicate potential technical problems.

CLSI recommends daily performance monitoring of fluorescence intensity, color compensation and verification of system performance using a QC specimen. Each laboratory must have QC procedures to monitor instrument performance. For hematology immunophenotyping, color compensation and verification of performance can be evaluated with each specimen panel.

Optics should not be changed by a clinical laboratory without specific validated procedures. It is more appropriate to have the manufacturer check and adjust optics during routine maintenance and to provide documentation for the laboratory^[3].

4.1 Fluorescence linearity Linearity of fluorescence detection should be checked on a monthly basis or as recommended by the instrument manufacturer. Linearity should be constant unless the laser is unstable. Monitoring linearity is a check of PMT voltages. Using the same PMT settings used for clinical specimens, a set mixture of four to five multi-level fluorescence beads, each with a known relative fluorescence intensity level, should be run. Products are available from several manufacturers, for example Flow Cytometry Standards, or Spherotech. Acceptable mean fluorescence intensity (MFI) ranges for each bead in the mixture should be established by 20 replicate runs over a 5 day time period. The correlation coefficient of the MFI versus fluorescent molecules per bead should be equal to or greater than 0.98 (1.0 is optimal). These same bead mixtures should then be run once a month and the linearity, or the relative difference in MFI between each of the beads, should remain constant for each fluorescence PMT. Graphs of linearity should be drawn following manufacturers' recommendations^[4,6].

4.2 Cross-instrument performance Laboratories that are performing the same clinical immunophenotyping protocols on more than one instrument should include a semi-yearly cross-instrument comparison of at least five different and representative clinical samples stained with each of the standard protocols used in the laboratory. The results obtained for each sample should not differ between instruments by more than a pre-defined acceptable variance. This cross-instrument sample testing process should be documented and corrective action plans must be established and followed when cross-comparison results fail to meet performance specifications.

5 Pre-analytical QC

Laboratories must not only establish acceptance criteria for immunophenotyping specimens, but must have procedures to assess acceptability and documentation to assist in troubleshooting and the interpretation of results. CLSI recommends a visual analysis on receipt to confirm specimen quality. Hemolysis, clot, partial draw (especially in ACD tubes), temperature extremes and improper labeling should be documented before transfer to the testing laboratory.

The laboratories must establish specimen requirements, recommended transport conditions and criteria for acceptability. All specimens must be considered infectious and must be handled in such a way as to minimize risks for the laboratory personnel. Procedures for handling, packaging, labeling and transporting potentially infectious biological specimens must be followed. Briefly, place specimen test tubes or vials (preferably unbreakable) in a leak-proof secondary container with sufficient absorbent material to absorb contents of the vessel, if breakage should occur. The U. S. -Canadian Consensus strongly recommended that every attempt be made to derive useful information from any specimen submitted for leukemia or lymphoma immu-

nophenotyping analysis. Requesting physicians should always be informed about sample quality issues; any compromised specimen, whether evaluated or not, must be described in the test interpretation.

Clinical specimens must be appropriately identified. The minimum information is a unique patient identifier (for cytogenetic analyses, two identifiers may be required), test ordered and date of sampling. Other information helpful for interpretation include presumptive diagnosis, age, gender, date and time of specimen collection, source of specimen, name of physician, and recent treatment or medications. Confidentiality must be assured and documentation tracking the specimen's handling is essential. Tests may only be performed if requested by an authorized person. Verbal requests must be followed by written authorization. All requests must be retained at least 2 years. The laboratory shall test leukemia/lymphoma specimens within 48 hours of collection^[4].

6 Reagent and method validation

An immunophenotyping procedure first requires the selection of monoclonal reagents and fluorochromes which are to be used in the multi-parameter analysis and interpretation of results. The hematopathology testing panel must be designed to clearly distinguish normal and abnormal immunoreactivity patterns based on differences in light scatter and/or fluorescence intensity. These pattern comparisons are integral to the interpretation of results. Therefore, the panel must include all relevant markers and reflect the visual expectations of the diagnostic interpreters.

The U. S. -Canadian Consensus Conference stopped short of recommending a diagnostic panel for leukemia and lymphoma immunophenotyping, however there was general agreement that at least 20 to 24 markers were required for complete characterization of acute leukemia. Similar numbers of markers were required for complete characterization of chronic lymphoproliferative disorders. Each specimen must be considered a unique case and fully evaluated to minimize missing any abnormality. The selection of panel reagents needs to balance the economic need for a minimal number of monoclonal antibodies to identify abnormal populations with the scientific need to detect abnormal antigen expression. The committee concluded that it is more important that each laboratory have adequate experience with the binding characteristics of their testing reagents than specifying what those reagents should be. How reagents are combined in testing panels is left up to the testing laboratory, but should be based on a target-oriented strategy^[7].

All the combinations have to be validated in each laboratory and the analysts must be familiar with the patterns associated with each combination. Different clones (under the same CD) may perform differently. Similarly, the same CD labeled with different fluorochromes will show a different level of intensity and possible color overlap. Fluorochrome intensity and antigen expression are both important in establishing the best antibody combination to evaluate antigen density. Some tandem dye combinations can help to minimize these problems.

Different combinations of antibodies can be used to provide the same clinical interpretation and no single panel will accommodate all leukemia or lymphoma phenotypes. To compound the difficulty, new markers and reagents, when determined to be clinically valuable, must be validated by each laboratory. Any panel of monoclonal reagents used by a laboratory must be chosen not only for their technical performance, but also to satisfy the experience and expectations of the technical and medical personnel who must differentiate abnormal populations from normal cells.

The proper combination of monoclonal antibodies within a cocktail must consider antigen expression on normal and abnormal cells as well as the fluorochrome combinations to minimize interference between reagents. For multiple fluorochrome combinations, the laboratory must not only validate that each monoclonal performs as expected in the procedure, but must also prove that the performance is not affected by co-blocking of epitopes, fluorescence quenching, or energy transfer.

Antibody cocktails may be prepared by the laboratory. Whether reagents are mixed immediately before each use or stored as a premix solution, performance of each cocktail must be validated and expiration dates documented.

Each monoclonal antibody in each cocktail must be titrated individually for optimal signal-to-noise separation. For each new multi-color combination added to the panel, it is necessary to first validate that the performance (mean intensity fluorescence and percent positive) of each antibody when used in combination is comparable (within 2 S. D.) to the performance of each antibody used alone at the same concentration on the same target cell population. If reagents are stored in a pre-mixed cocktail format, it is also important to institute sufficient QC procedures, using appropriate control cells, to verify stability of performance of the cocktailed reagents over time. Stable cocktails will provide results wherein the mean channel and percent of gated cells will not differ more than 2 S. D. between cocktailed and single reagents. Documentation must be summarized and maintained in the laboratory and should include instrument PMT settings and compensation. No more than one combination of PMT and compensation may be used within a panel. The best approach is to set the PMT with the single reagent and then check compensation with the final cocktail^[5-6].

7 Sample preparation

Sample preparation for flow cytometric analysis must consider the type of specimen submitted and the number of cells available for analysis. Peripheral blood, bone marrow, or fresh tissue specimens should be processed to maintain all pertinent cellular and antigenic parameters for lymphoid, monocytic and myeloid cells, eliminating erythrocytes, at a concentration optimal for monoclonal staining ($0.2^{-2} \times 10^7/\text{mL}$; 0.1 mL per tube)^[3]. All processing procedures must be written, approved and daily records maintained. Test records must be retained or at least 2 years.

Integral to the analytical procedure, the laboratory must decide on and validate a sample preparation procedure. Erythro-

cyte lysis is the most commonly used method for preparing peripheral blood or bone marrow specimens for clinical immunophenotyping. Lysis can be performed either before staining with monoclonal antibodies or after. Whichever method for erythrocyte lysis is used must be validated prior to reporting clinical results. Lysis is optically monitored by observing a cloudy red suspension of cells change into a clear, translucent solution in 5 to 10 min. If red cells are still present after centrifugation and resuspension, the lysing process should be repeated.

Tissue specimens should be disaggregated and filtered to create a single cell suspension prior to staining. This laboratory procedure must also be written and validated. Cytological slides should be prepared and reviewed with the H&E slides of the tissue to assess the loss of fragile cells, such as Reed-Sternberg cells. Some large cells in fibrotic tissues may be difficult to remove and cell loss will occur during the processing. Other sampling errors are due to focal lymphoma in the lymph node, not seen in the isolated cells or in the fixed tissue. Cell fragility may be minimized by gentle tissue processing and teasing the tissue apart with a needle.

Reproducible monoclonal antibody staining relies on the proper proportion of antibodies to cells. For hematolymphoid immunophenotyping, cell counts should be performed on all specimens to ensure correct proportions of cells to monoclonal antibodies and thus accurate staining intensity patterns. CLSI recommends that normal WBC counts in the sample should range from 0.2 to 2×10^6 nucleated cells per tube. (There are no similar recommendations for bone marrow and each laboratory should establish and validate its own procedures.) Light scatter patterns must be appropriate after lysing and staining^[3].

7.1 Viability Assessment of viability is crucial for evaluating leukemia and lymphoma specimens by flow cytometry because cell membrane integrity is required for antigen expression. Viability at the time of collection, the influence of transportation environment and the time before testing can all affect sample viability and compromise tumor cell detection^[3-4].

Flow cytometry performed on leukemias and lymphomas should be held to less stringent viability restrictions than quantitative immunophenotyping. When viability is low, for example less than 50%, a gating discrimination strategy should be used to limit the analysis to live cells, minimizing nonspecific binding from dead cells which could lead to misdiagnoses. Several different methodologies have been used to report viability, trypan blue, 7-aminoactinomycin D (7-AAD), or propidium iodide (PI). The most common is trypan blue dye exclusion visualized microscopically with a hemocytometer. Dead cells stain blue because their membranes are broken, allowing the dye to enter the cell. Viability is reported as a percentage of the total cell population.

The exclusion of dead cells is now being commonly performed on the flow cytometer using dyes such as 7-AAD. Dead cell discrimination can be performed simultaneously with im-

munophenotyping. The most common three-color application in single laser instruments is FITC in FL1, PE in FL2 and 7-AAD in FL3. In this setting, 7-AAD (Emax = 655 nm) is preferred over PI (Emax = 625 nm) because it has less spectral overlap with PE (Emax = 578 nm).

7.2 Staining process QC Once the instrument setup and performance validation is complete, a normal blood is run to check the staining process^[4]. New York State recommends that a freshly prepared whole blood specimen from a healthy donor be included as a positive control; at least monthly for hematolymphoid immunophenotyping. For this purpose, there are some commercially available preserved cell preparations, such as CD-Chex, CD-Chex Plus (Streck Labs), Immunotrol (Beckman Coulter), which contain most of the cellular markers evaluated in hematolymphoid immunophenotyping. Staining each new lot of stabilized cells in duplicates will provide a range to be used for quality control. Most importantly, each leukemia or lymphoma specimen will contain a normal cell population which can be used to document appropriate staining patterns.

7.3 Isotype controls The appropriate use of isotype controls is still controversial in the clinical flow cytometry community. For quantitative immunophenotyping, the isotype control was included to provide a negative cell reference and to set up the axes for assessing the numbers of cells in each quadrant. However, immunophenotyping of leukemias and lymphomas is a qualitative assessment and does not usually require quadrant analysis. Furthermore, since all leukemias and lymphomas contain a mixture of cells, those cells that do not express the antigen can be served as the negative control for the ones that do.

Most of the commercially available antibodies have been selected by screening for the IgG class of immunoglobulin. For IgG class antibodies, however, it is important to recognize that some IgG subclasses will be more problematic due to their increased binding to Fc receptors present on various cell types. In general, the order of 'stickiness' is IgG2b > IgG2a > IgG1. As a result of their carbohydrate structure, some antigenic targets (e.g., CD15) will invariably result in generation of only IgM class antibodies. High non-specific binding of IgM class antibodies is commonly observed and can be influenced by the donor, the target cell type, and the particular staining method used. For cytoplasmic staining, the use of isotype controls is recommended due to the significant influence of cell size on the degree of non-specific staining. It is important that each monoclonal antibody has the corresponding fluorochrome and isotype reagent to monitor non-specific binding performance.

The presence of non-specific binding of an isotype can affect the interpretation of the fluorescence intensity patterns and, as such, is critically important in the analysis of hematolymphoid immunophenotyping.

7.4 Procedure control Fluorescence staining intensity patterns are more easily interpreted if there is a common factor in each tube. For example, CD45 can be used as the 'anchor' in each tube of a multi-color cocktail, providing a reference population for comparisons to other antigenic expression patterns.

The most important parameter, however, is experience-based. The analyst must understand the markers and the patterns used to identify an abnormal population. Running a normal control will provide a baseline for all the markers. Any deviation from the normal pattern should be reported and confirmed with other markers or technologies, such as cytogenetic or molecular techniques.

8 Data analysis and interpretation

Hematolymphoid immunophenotyping by flow cytometry can be used to distinguish abnormal cell populations from normal cells and, in the process, identify the profile of antigens in the abnormal population. Various patterns of immunoreactivity with monoclonal antibodies are combined with light scatter to characterize populations of interest. The specific reagent combinations can be chosen to identify abnormal antigen expression and to characterize cells at various stages of maturation.

The complement of monoclonal antibodies necessary to identify hematolymphoid malignancies is not standardized^[1,7]. In fact, strategies for immunophenotyping must consider antibody specificity, fluorochrome, laboratory workload and, most importantly, appropriate multi-parameter combinations. Most laboratories use a comprehensive screening panel and add confirmatory monoclonal antibodies if needed. Reimbursement for laboratory testing is a major factor, thus laboratories strive to make diagnoses with a minimal number of antibodies. The U. S. -Canadian Consensus Conference estimated that the average North American laboratory uses 19 antibodies to diagnose leukemia and 16 for lymphoma.

Gating is one of the most important parameters in multi-parameter data analysis. In a diagnostic case, all cells should be evaluated. The U. S. -Canadian consensus recommended that the initial analysis include all viable cells. Further analyses on gated populations should only be performed if all cells of interest are contained within the analysis gate. Strategies for specific diseases may include gating on B cells to determine clonality or identification of leukemic blasts in a CD45 versus SSC gating strategy^[7-8].

Data analysis of abnormal populations requires that fluorescence intensity be measured. The calculation of percent positive cells does not aid in the interpretation, while presence of inappropriate determinants or inappropriate intensity expression of antigens can be diagnostic. Resolving dimly positive populations from negative populations is crucial in the accurate assessment of fluorescence intensity. There are no hard and fast rules for defining these populations; in fact, reagent combinations are generally optimized to provide discrimination between these populations based on the subjective review of the pathologist. For critical markers, such as immunoglobulin light chains, in the diagnosis of non-Hodgkin lymphoma, the redundant use of antibodies to light chains using alternative fluorochromes, anti-

body pairs from different vendors, or monoclonal versus polyclonal reagents, can significantly enhance confidence in interpretation of monoclonality in the face of very low level expression^[5-7].

The interpretation of the flow cytometric data relies on an experienced diagnostic interpreter, usually a hematopathologist, reviewing all data, including morphology. Drawing conclusions from either microscopy or flow cytometry is an experience-based skill and the final interpretation should reflect a synthesis of all available information about the case. All reports must be sent to the requesting party and must also be retained for a minimum of 2 years^[9-10].

参考文献

- [1] Jaffe ES, Harris NL, Stein H, et al. World Health Organization Classification of Tumors, Pathology and Genetics of Tumors of Haematopoietic and Lymphoid Tissues[M]. IARC Press: Lyon, 2001.
- [2] Swerdlow SH, Campo E, Harris NL, et al. WHO Classification of Haematopoietic and Lymphoid Tissues[M]. IARC Press: Lyon, 2008.
- [3] Stelzer GT, Marti G, Hurley A, et al. U. S. -Canadian Consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry; standardization and validation of laboratory procedures[J]. Cytometry, 1997, 30(5): 214-230.
- [4] Borowitz MJ, Bray R, Gascoyne R, et al. U. S. -Canadian Consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry; data analysis and interpretation[J]. Cytometry, 1997, 30(5): 236-244.
- [5] Braylan RC, Atwater SK, Diamond L, et al. U. S. -Canadian Consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry; data reporting[J]. Cytometry, 1997, 30(5): 245-248.
- [6] Davis BH, Foucar K, Szczarkowski W, et al. U. S. -Canadian Consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry; medical indications[J]. Cytometry, 1997, 30(5): 249-263.
- [7] Hurley AA. Quality control in phenotypic analysis by flow cytometry[J/OL]. Curr Protoc Cytom, 2001-05-01 [2012-08-26], <http://www.ncbi.nlm.nih.gov/pubmed/18770714>.
- [8] Clinical and Laboratory Standards Institute. H43-A2 Clinical cytometric analysis of neoplastic hematolymphoid cells; approved guideline-second edition[S]. Clinical and Laboratory Standards Institute, Pa, USA: 2007.
- [9] Carey JL, McCoy JrP, Keren DF. Flow cytometry in clinical diagnosis[M]. 2nd Edition 2008 American Society for Clinical Pathology, Press, Chicago IL.
- [10] New York State Department of Health. Clinical Laboratory Standards of Practice: Cellular Immunology[S/OL]. 2012-06-25 [2012-08-13]. <http://www.wadsworth.org/labcert/lep/standards.htm>.